

Health AI- Intelligent Healthcare

Assistant

Generative AI with IBM



1.INTRODUCTION:

Health AI – Intelligent Healthcare Assistant

Team Members

- 1. K. Mahendrabupathi - (88472D5E26B209EA06B4C6E173DE899B)**
- 2. R. Pradeepraja - (AC97844F35595E9813B44B55BB6D70A9)**
- 3. P. Murugesan - (E3D4E9CF22BF393834542FB2C1F30F77)**
- 4. S. Naresh - (3C5B5E72979BBE6B9A04C946801155D5)**

2.PROJECT OVERVIEW :

Health AI-intelligent healthcare Assistant is an AI-driven system designed to improve healthcare delivery through intelligent support. The solution integrates AI/ML models with LLMs and vector databases, user-friendly healthcare assistance.

• Conversation Interface

The conversational interface allows users to interact with the healthcare natural text or voice. It supports health queries, policy guidance and human experts.

• Policy summarization

The system automatically extracts and condenses healthcare policies into simple,

easy-to-read summaries. It helps patients, citizens and staff quickly understand key update.

- **Eco-Tip Generator**

The eco-tip generator provides simple, actionable health and environment friendly suggestions. It promotes sustainable practices like waste reduction, energy, saving methods for community health.

- **Citizen feedback loop**

The citizen feedback loop collects patient and public opinion through survey, chat, or forms. It analyses feedback using sentiment and topic detection to improve system enhancement.

- **KPI Forecasting**

KPI forecasting predicts key healthcare performance indicators such as bed occupancy, patients' inflow, staff availability, and resource usage. It helps in planning, decision-making, improving overall Healthcare efficiency.

- **Anomaly Detection**

Anomaly detection identifies unusual pattern in healthcare data, such as sudden diseases spikes or abnormal resources usage. It enables early alerts and quick action to improve patient safety.

- **Multimodal input support**

Multimodal input support allows a healthcare assistant to receive and process text, voice, image and sensor data. This enables more natural, accurate, and efficient interaction for patients and medical staff.

- **Stream lit to Gradio UI**

Stream lit to Gradio UI allows converting python apps into interactive web interface with easy deployment. Gradio provide simple drag-and-drop components for inputs and outputs, user interaction.

3.ARCHITECTURE:

- Front-end Architecture of Health AI-Intelligent Healthcare assistant interaction such as web, mobile or voice app, allowing patients or citizens to accessible way.
- Back-end Architecture of Health AI-Intelligent Healthcare assistant handles business logic, connect LLM for prediction and recommendations provide APLs for secure, real-time communication with the front end
- **LLM Integration:**

The LLM Processes user queries, interpret intent, and generate response while interacting with for knowledge retrieval and personalized recommendations.

- **Vector sector:**

The vector sector in health care AI stores embedding of medical data, policies, records to enable fast semantic search, allowing the AI to Retrieve contextually relevant information guidance.

- **ML Modules:**

The ML Modules provide analytics, anomaly detection and personalized health recommendations by analysing patient data, trends medical knowledge to support.

4.SETUP INTRODUCTION:

- **Prerequisites:**

- Install python, node and code editor
- Access open AI API and vector database
- Basic ML knowledge and libraries
- Front end skills (react, angular or flutter)
- Medical database or pre-trained models

- **Installation process:**

- Install python, node and set up a virtual environment.
- Install backend and AI libraries
- Setup vector database and create index.
- Run back end and front end test queries and AI responses.

5.FOLDER STRUCTURE:

- App: application code
- Data: Storage data files
- Docs: Documentation
- Test: Unit test
- Ven: Virtual environment
- Txt: Text
- Utils: Utility function

6.RUNNING THE APPLICATION:

To start the project:

- Open terminal and clone the project from GitHub.
- Navigate to the project folder

- Health AI-Intelligent healthcare assistant-AI.
- Create and activate a virtual environment.
- Install dependencies using pip install -r requirements. Txt.
- Set up database and vector store (initialize configs).
- Configure API keys (OpenAI/LLM, database, etc.).
- Start backend server using python backend/ API/main.py.
- Run frontend UI using stream lit run frontend/ app. Py Or gradio app.

Py. [Frontend \(Stream lit\):](#)

Enter frontend folder - run stream lit – open browser- use Health AI-Intelligent healthcare assistant UI.

[Backend \(fast API\):](#)

Enter backend folder – run fast API server – open API docs – backend ready for health AI-Intelligent healthcare assistant UI integration.

7.API DOCUMENTATION:

API Documentation provides endpoint for patient interaction, health record management, and AI – drive recommendations, enabling seamless healthcare support through conversational and analytical modules.

8.AUTHENTICATION:

Register user/client – user signs up or application register to receive API credentials. Generate API key/token – system issues a unique API key or JWT token. Send authentication request – client sends login request with credentials. Receive access token – server validates and responds with an access token. Use token in API calls – includes authorization: bearer <token> in header for every request.

9.USER INTERFACE:

- Provides a simple login for secure access.
- Chat window support text and voice queries.
- Dashboard shows health records and reports.
- AI gives personalized tips and recommendations.

10.TESTING:

- Testing in the health AI – intelligent healthcare assistant ensures the system works reliably and safety for patients and doctors. First, unit testing is

performed to verify individual modules like authentication, chat, and health records access.

- User acceptance testing (UAT) then validate that real user finds the system useful, accurate, easy to use. This step-by-step testing process guarantees a robust and trustworthy healthcare assistant.

11.KNOWN ISSUES:

The Health AI-Intelligent healthcare assistant, while highly efficient, currently faces certain limitations. Occasional inaccuracies in symptom analysis may occur due to limited training data for rare conditions. Multimodal input processing, such as interpreting images and text together, may experience delay or misinterpretation under poor-quality inputs. Some users may encounter minor UI inconsistencies across different devices. These issues, improve reliability and enhance overall user experience. During testing and initial deployment, several issues were identified that may affect the system performance or user experience. Some features may respond slowly under heavy load and occasionally UI observed on different devices. Certain data inputs may trigger unexpected errors if they do not match the expected format. Additionally, integration with external APIs can cause server downtime, which will be addressed in future releases to ensure stability and functionality.

12.FUTURE ENHANCEMENTS:

In the future, as the Health AI-Intelligent healthcare assistant updates the system, improvements for the health AI assistant will focus on expanding its capabilities and improving patient care. Planned enhancements include incorporating more advanced diagnosis tools, supporting multimodal input such as voice and image, and integrating with additional health data sources for personalized recommendations. The system will also implement predictive analytics for early disease detection, real-time health monitoring and enhance natural language understanding for more accurate and empathetic interactions. The user interface will be refined for greater accessibility and responsibility while automated error handling and notification systems will be implemented to ensure smoother operations. These enhancements aim to provide a more robust, efficient and user-friendly solution.

13.PROJECT SCREENSHOT:

```
import gradio as gr
import torch
from transformers import AutoTokenizer, AutoModelForCausalLM

# Load model and tokenizer
model_name = "The-Grazer/gradio-1.2.0-instruct"
tokenizer = AutoTokenizer.from_pretrained(model_name)
model = AutoModelForCausalLM.from_pretrained(
    model_name,
    torch_dtype=torch.float16 if torch.cuda.is_available() else torch.float32,
    device_map="auto" if torch.cuda.is_available() else None
)

if tokenizer.pad_token is None:
    tokenizer.pad_token = tokenizer.eos_token

def generate_response(prompt, max_length=512):
    inputs = tokenizer(prompt, return_tensors="pt", truncation=True, max_length=512)

    if torch.cuda.is_available():
        inputs = {k: v.to(model.device) for k, v in inputs.items()}

    with torch.no_grad():
        outputs = model.generate(
            **inputs,
            max_length=max_length,
```

```
max_length=max_length,
temperature=1,
do_sample=True,
pad_token_id=tokenizer.eos_token_id
)

response = tokenizer.decode(outputs[0], skip_special_tokens=True)
response = response.replace(prompt, "").strip()
return response

def disease_prediction(symptoms):
    prompt = f"Based on the following symptoms, provide possible medical conditions and general medication suggestions. Always emphasize the importance of consulting a doctor for personalized advice."
    return generate_response(prompt, max_length=512)

def treatment_recommendation(age, gender, medical_history):
    prompt = f"Generate personalized treatment suggestions for the following patient information. Include how medication and general medication guidelines/medical conditions."
    return generate_response(prompt, max_length=512)

# Create simple interface
with gr.Blocks() as app:
    gr.Markdown("Medical AI Assistant")
    gr.Markdown("Disclaimer: This is for informational purposes only. Always consult healthcare professionals for medical advice.")

    with gr.Tab():
        with gr.Tab("Disease Prediction"):
            with gr.Text():
                symptoms_input = gr.Textbox()
```

```
symptoms_input = gr.Textbox(
    label="Enter symptoms",
    placeholder="e.g., fever, headache, cough, fatigue...",
    lines=1
)
predict_btn = gr.Button("Predict Disease")

with gr.Column():
    prediction_output = gr.Textbox(label="Possible conditions & recommendations", lines=10)

predict_btn.click(disease_prediction, inputs=symptoms_input, outputs=prediction_output)

with gr.Tab():
    with gr.Tab("Treatment Plan"):
        condition_input = gr.Textbox(
            label="Medical condition",
            placeholder="e.g., diabetes, hypertension, migraine...",
            lines=1
        )
        age_input = gr.Number(label="Age", value=30)
        gender_input = gr.Dropdown(
            choices=["Male", "Female", "Other"],
            label="Gender",
            value="Male"
        )
        history_input = gr.Textbox(
            label="Medical history",
```


Medical AI Assistant

Disclaimer: This is for informational purposes only. Always consult healthcare professionals for medical advice.

[Disease Prediction](#) [Treatment Plans](#)

Enter Symptoms

None

Analyze Symptoms

Possible Conditions & Recommendations

1.6. Immune-related or autoimmune disorders: Conditions like lupus, rheumatoid arthritis, or other autoimmune diseases may present with fever, in addition to other systemic symptoms.

General medication suggestions:

2. Acetaminophen or Ibuprofen: These pain relievers and fever reducers can help alleviate fever and associated discomfort, however, always consult your doctor before starting any medication, especially if you have pre-existing medical conditions or are taking other medications.

3. Antibiotics: If a bacterial infection is suspected based on symptoms and examination, antibiotics might be prescribed.

4. Antiviral medications: For specific viral infections like influenza or COVID-19, antivirals may be recommended by a healthcare professional.

5. Hydration and rest: Taking extra fluids and ensuring adequate sleep can support overall recovery.

6. Monitoring symptoms: Keep track of the severity and duration of symptoms, as well as any worsening of underlying conditions or new symptoms, to assist in identifying potential complications.

In summary, while fever can be a symptom of various conditions, the most common causes involve viral infections. To effectively manage symptoms, always consult a healthcare professional for a proper diagnosis and personalized treatment plan.

View in PDF | Build with Gemini | Settings

Medical AI Assistant

Disclaimer: This is for informational purposes only. Always consult healthcare professionals for medical advice.

[Disease Prediction](#) [Treatment Plans](#)

Enter Symptoms

None

Analyze Symptoms

Possible Conditions & Recommendations

General medication suggestions:

2. Acetaminophen or Ibuprofen: These pain relievers and fever reducers can help alleviate fever and associated discomfort, however, always consult your doctor before starting any medication, especially if you have pre-existing medical conditions or are taking other medications.

3. Antibiotics: If a bacterial infection is suspected based on symptoms and examination, antibiotics might be prescribed.

4. Antiviral medications: For specific viral infections like influenza or COVID-19, antivirals may be recommended by a healthcare professional.

5. Hydration and rest: Taking extra fluids and ensuring adequate sleep can support overall recovery.

6. Monitoring symptoms: Keep track of the severity and duration of symptoms, as well as any worsening of underlying conditions or new symptoms, to assist in identifying potential complications.

In summary, while fever can be a symptom of various conditions, the most common causes involve viral infections. To effectively manage symptoms, always consult a healthcare professional for a proper diagnosis and personalized treatment plan.

View in PDF | Build with Gemini | Settings

Medical AI Assistant

Disclaimer: This is for informational purposes only. Always consult healthcare professionals for medical advice.

[Disease Prediction](#) [Treatment Plans](#)

Medical Condition

diabetes

Age

41

Gender

Female

Medical History

allergies

Generate Treatment Plan

Personalized Treatment Plan

1. Dietary Modifications:

- Include more whole foods like fiber, lean proteins, and healthy fats like fruits, vegetables, legumes, nuts, and seeds. Examples include oatmeal with berries and nuts for breakfast, grilled chicken salad with olive oil dressing, and bowls of Greek yogurt with honey and almonds.
- Limit intake of refined carbohydrates, sugary snacks, and sugary beverages like sodas. Opt for whole grains instead of refined grains.

2. Regular Exercise:

- Aim for at least 150 minutes of moderate intensity or 75 minutes of high-intensity aerobic activity per week, along with strength training exercises on 2 or more days a week.
- Consider low-impact activities like swimming or cycling if joint stress is a concern.


3. Weight Management (if applicable):

- If overweight, consider losing 5-10% of body weight through a combination of dietary restrictions and increased physical activity.
- Aim for a healthy BMI (Body Mass Index) of 18.5-24.9.

4. Monitor Blood Glucose Levels:

- Regularly check blood glucose levels using a home glucose monitoring system, following your healthcare provider's guidelines.
- Keep a record of readings, especially before and after meals, to identify patterns and adjust treatment as needed.

View in PDF | Build with Gemini | Settings



14.VIDEO LINK:

<https://drive.google.com/file/d/1tfZugVNZsBI8KKwDaEexnMUk5M2wpYI/view?usp=sharing>

THANK YOU